

HEALTH IS THE HEART OF OUR BUSINESS.

HOLTER-RECORDER EP830 / EP830-12

DIGITAL HOLTER RECORDER WITH PACEMAKER DETECTION



High patient acceptance

3 or 12 channel recording

Long term use with more than 14 days of continuous 3-channel recording

OLED color display showing 3 channel real time ECG simultaneously

Protection against splash water & dust



MANAGE COMPLEX DATA EASILY

AMEDTEC ECG^{pro}

TECHNICAL DATA

Dimensions	Size	80 x 59 x 18 mm (L x W x H)
	Weight	63 g
	Operating position	any orientation
Electrical	Patient cable	4, 5, or 10 wires
	Battery type	1 AA battery: Alkaline, Lithium or NiMH
	Applied part	Type CF
Functional	Recording channels	3 or 12
	Sensitivity	2.5 μ V per LSB
	AD resolution	24 bit
	Recording	full disclosure
	Data transfer	USB 2.0 or card reader
	Sample rate	250/s, 500/s or 1000/s
	Frequency response	0.05 to 240 Hz
	Display	OLED color display
	PM detection	permanently
	Waterproof	IPX6 grade
	Recording time	3 channel recording: > 14 days 12 channel recording: 2 days


ACCESSORIES



Holter ECG electrodes



disposable pouch




disposable sanitary
neck strap



patient cable

AMEDTEC Medizintechnik Aue GmbH
Schneeberger Straße 5 · 08280 Aue
Germany
Phone: +49 (0)3771 59 82 70
Fax: +49 (0)3771 59 82 790
info@amedtec.de · www.amedtec.de



CE Certificate:
0494 DIN EN ISO 13485

Errors and omissions excepted
M0703 019-1942-2_01